



William G. Kaylakie DMD
11615 Angus Road, Suite 113
Austin, TX 78759
(512) 794-2822

Orthotic Splint

- Brush splint with denture tablet dissolved in water, or with mild detergent and water. **Do not use toothpaste.**
- To clean, place splint in a cup with denture tablet or powder and water. After 15 to 20 minutes, rinse the splint well.
- Splint can be stored in a cup or retainer case and a small amount of mouthwash, if desired. **(DO NOT store splint in retainer case submersed in water.)**